Instructions for the use of IGS Instant Grip Socks

IGS Instant Grip Socks are easy to use, breathable and flexible. The socks prevent slipping while walking on floors and wet surfaces and ease changing position in bed. The product is made of environmentally friendly material that can be recycled for energy. The Instant Grip Socks create strong friction between the patient's feet and the floor. When dressing the patient, roll the sock up foot-do not pull.

Selection of IGS Instant Grip Sock

Choose the appropriate sock size (S, M, L, XL, XXL) based on the user's shoe size. If there are swelling, sutures or ulcers on the foot, simply go up one size. You can stretch the sock out at points of swelling, sutures or ulcers as required by cutting holes in the mesh structure. The mesh will not come apart beyond the cuts made. Not recommended for patients with foot drop or other neurological issue affecting gait.

Dressing

1. Start by placing the tip of the rolled sock at the toes and continue by unwrapping the sock towards the heel- do not pull.
2. Roll the sock open towards and over the ankle by gently stretching the sock with your thumbs.
3. Level the sock's mesh structure by "rubbing" it on the foot, at the same time checking that the sock is not too tight at the toes.

NB! Avoid straining the mesh structure at the heel and avoid tight layers at the ball of the foot and ankle.

Cutting

4. If the user's feet are swollen or there are ulcers or stitches, the sock can be widened by cutting. Raise the mesh off the skin and cut at the chosen locations.
5. To wash toes in the shower, you can cut approximately 5cm of the sock horizontally at the lower part of the ball of the foot. Flip the tip of the sock under the foot and after washing turn it back into place.

NB! Always leave at least 2cm of intact material at the sock stem in order to keep the sock in place.

Undressing

6. Remove the socks by rolling from the cuff of the sock down towards the tip.

After using the Instant Grip Socks for 24 hours on bare feet the socks need to be rinsed or washed in 40 degrees with mild detergent.

The product is not suitable for outdoor use.

Not recommended for people with strong rubber allergies or on bare ulcerated skin.

Not recommended for people who drag their legs

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CareCare Oy