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
care  
care  
www.carecare.fi

## Using the Grip Sock

The Grip Sock is easy to use, flexible and breathable. Prevents slipping on dry and damp surfaces and makes it easy to change position in bed. Made of environmentally friendly, recyclable material. Put it on by rolling, not pulling from the stem. Do not use outdoors.

### Sock selection

Choose a sock according to the size of the user's foot (XXS-XXL). For a leg with swelling, stitches, or ulcers, choose a larger size. If necessary, you can loosen the sock by cutting the openings in the mesh structure larger without unravelling. Protect ulcerated skin with wound dressings before using the sock. The sock can penetrate substances that increase the slipperiness of the floor surface (e.g. silicone in a wound dressing). Not recommended for people with severe rubber allergies or foot lagging.

### Dressing

1. Start by placing the tip of the rolled sock at the toes and continue by unwrapping the sock towards the heel- do not pull.
2. Roll the sock open towards and over the ankle by gently stretching the sock with your thumbs.
3. Level the sock's mesh structure by "rubbing" it on the foot, at the same time checking that the sock is not too tight at the toes.

NB! Avoid straining the mesh structure at the heel and avoid tight layers at the ball of the foot and ankle.

### Cutting

4. If the user's feet are swollen or there are ulcers or stitches, the sock can be widened by cutting. Raise the mesh off the skin and cut at the chosen locations.
5. To wash toes in the shower, you can cut approximately 5 cm of the sock horizontally at the lower part of the ball of the foot. Flip the tip of the sock under the foot and after washing turn it back into place.

NB! Always leave at least 2 cm of intact material at the sock stem in order to keep the sock in place.

### Undressing

6. Remove the socks by rolling from the cuff of the sock down towards the tip.

After using the Instant Grip Socks for 24 hours on bare feet the socks need to be rinsed or washed in 40 degrees with mild detergent.

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